

# X-Gal Staining™

amsbio

## Assay Kit

| Cat. No.                               | Content   | Qty      | Storage |
|--|---|----------|---------|
| A10300K<br>(50 assays in 60 mm dishes) | Fixing Buffer   | 125 ml   | 4°C     |
|  | 10X PBS   | 75 ml    | 4°C     |
|  | Staining Buffer   | 125 ml   | 4°C     |
|  | 25X X-Gal Stock (5-bromo-3-indolyl- $\beta$ -D-galactopyranoside) | 4 X 1 ml | -20°C   |

|                               |  |
|-------------------------------|--|
| <b>Shipping &amp; Storage</b> | Shipped on blue ice. The X-Gal Stock solution should be stored at -20°C. All other components should be stored at 4°C. Stable for 6 months when stored properly. |
|-------------------------------|--|

| RELATED PRODUCTS   | Catalog # |
|--|-----------|
| GenePORTER™ Transfection Reagent, 75 reactions             | T201007   |
| GenePORTER™ Transfection Reagent, 150 reactions            | T201015   |
| GenePORTER™ Transfection Reagent, 750 reactions            | T201075   |
| GenePORTER™ 2 Transfection Reagent, 75 reactions           | T202007   |
| GenePORTER™ 2 Transfection Reagent, 150 reactions          | T202015   |
| GenePORTER™ 2 Transfection Reagent, 750 reactions          | T202075   |
| gWiz™ $\beta$ -galactosidase Expression Vector, 25 $\mu$ g | P010200   |
| Enhanced $\beta$ -galactosidase Assay Kit (CPRG)           | A10100K   |
| $\beta$ -Galactosidase Staining Kit (ONPG)                 | A10200K   |

**INTRODUCTION** *LacZ* is a commonly used reporter gene in transfection experiments because the gene product,  $\beta$ -galactosidase ( $\beta$ -gal), is very stable and resistant to proteolytic degradation and easily assayed. This assay kit provides all the required reagents, and offers a rapid and simple method to determine the percentage of cells transfected with *LacZ* expressing plasmids, such as Genlantis' gWiz  $\beta$ -gal vector.  $\beta$ -gal catalyzes the hydrolysis of  $\beta$ -galactosides (i.e. X-Gal) and cells transfected with  $\beta$ -gal expressing plasmid appear blue following fixation and incubation with X-Gal substrate. Blue cells can be visualized by microscopy.

- USAGE:**
- Transfect cells with a plasmid expressing *LacZ* gene.
  - Fix the cells with formaldehyde-glutaraldehyde buffer.
  - Stain the cells with X-Gal staining solution.
  - Observe the cells with blue stain under a microscope.
  - Calculate the percentage of stained cells in the total population versus non-transfected cells

## EXPERIMENTAL PROTOCOL

### A. Buffer Preparation

1. Dilute 10X PBS to 1X with distilled deionized water before use. 1X PBS may be stored at 4°C or room temperature for future use.
2. Dilute 25X X-Gal stock to 1X with Staining Buffer. Discard unused 1X X-Gal.

### B. Assay Protocol

Use the following table for recommended buffer volumes to use depending on the type and size of your tissue culture plate or dish:

| Type of culture dish | Fixing Buffer ( $\mu$ l/well) | Staining Buffer ( $\mu$ l/well) | 1X PBS Washing Buffer ( $\mu$ l/well/wash) |
|----------------------|-------------------------------|---------------------------------|--|
| Chambered slide      | 500                           | 500                             | 1000                                       |
| 24-well plate        | 250                           | 250                             | 500  |
| 12-well plate        | 500                           | 500                             | 1000                                       |
| 6-well plate         | 1000                          | 1000                            | 2000                                       |
| 60 mm dish           | 2500                          | 2500                            | 3000                                       |
| 100 mm dish          | 5000                          | 5000                            | 8000                                       |

1. Aspirate the medium 24-72 hours after transfection from the culture dish.
2. Wash the cells 1 time with 1X PBS.

3. Add Fixing Buffer to the dish and incubate for 10-15 minutes at room temperature.

**CAUTION:** Fixing Buffer contains chemicals that are corrosive, carcinogenic, and toxic. Handle Buffer carefully (see Materials Safety Data Sheet for further details) by wearing gloves, goggles, lab coats, and protective gear.

4. Remove the fixing solution from the dish and gently wash the cells 2 times with 1X PBS.
5. Add freshly prepared 1X X-Gal staining solution to the dish. Incubate the cells between 1-18 hours at 37°C in a humidified incubator. Adjust the incubation time according to the transfection efficiency.
6. Remove the X-Gal staining solution and wash the cells 1 time with 1X PBS.
7. Add 1X PBS to the dish. Examine the dish under a light microscope; count the stained and unstained cells in randomly selected fields. Calculate the percentage of stained cells in the total population.
8. To store the plates for weeks or months, fix each well with 1ml of 10% formalin in PBS (not supplied) for 10 minutes at room temperature. Rinse with 1X PBS and store in 1X PBS or 70% glycerol solution (not supplied) at 4°C.



#### UK & Rest of World

184 Milton Park, Abingdon  
OX14 4SE, Oxon, UK  
Tel: +44 (0) 1235 828 200  
Fax: +44 (0) 1235 820 482

#### Switzerland

Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
Tel: +41 (0) 91 604 55 22  
Fax: +41 (0) 91 605 17 85

#### Deutschland

Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
Tel: +49 (0) 69 779099  
Fax: +49 (0) 69 13376880

#### North America

23591 El Toro Rd, Suite #180  
Lake Forest, CA 92630  
Tel: +1 800 987 0985  
Fax: +1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com  
AMS Biotechnology