

## BG01V

Human embryonic stem cells.

<b>Catalog No.</b>	GSC-1103
<b>Cell Source</b>	Human – karyotypic variant (BG01V)
<b>Product Unit</b>	Cryopreserved vial
<b>Size</b>	~2 million cells/vial
<b>Doubling Time</b>	~36 hours
<b>Medium Renewal</b>	Every day
<b>Complete Growth Medium</b>	20% Serum Replacer in ES-DMEM/F12 (GSM-1002) plus bFGF (GSR-2001)
<b>Subculture Frequency</b>	Every 4-5 days
<b>Cryopreservation Medium</b>	hESfreeze (GSM-4200)
<b>Storage</b>	Vapor phase of liquid nitrogen
<b>Shipping</b>	Dry ice
<b>Safety Information</b>	<p><b>BioSafety Level: 1</b></p> <p>Safety conditions pertaining to well characterized agents not known to cause disease in healthy adult humans and of minimal potential hazard to laboratory personnel and the environment should be used. Lab coats and gloves should be worn at all times. Face protection is recommended when thawing a frozen vial stored in a liquid nitrogen freezer.</p>

### Intended Use

BG01V cells are derived from the wild-type BG01 human embryonic stem cells and are highly characterized. These cells are trisomy for chromosomes 12 and 17 and are XXY. Colonies of BG01V grown on irradiated or mitomycin C-treated mouse embryonic fibroblasts (MEFs) exhibit a predictable growth rate, form large, uniform colonies and are pluripotent. BG01V can be grown with feeders or feeder-free in the proper culture conditions. For Research Use Only (RUP). **Caution:** Not intended for human or animal diagnostic or therapeutic uses.

### Directions for Use

*All medium and reagents used in the culture of this product should be warmed to 25–37°C before use. Perform all activities under aseptic culture conditions.*

1. Prepare tissue culture dish with coating or feeder layer prior to thawing cells.
2. Place the frozen vial into a 37°C water bath as soon as possible and retrieve the vial before the contents are completely thawed (1–2 minutes).
3. Immediately transfer the contents of the vial to a 15-mL tube and dilute 1:10 with complete growth medium.
4. Spin the tube at 270 x g for 5 minutes in order to pellet the cells.
5. Resuspend the pellet in complete growth medium and perform a cell count.
6. Plate the cells to a T-25 flask or 6 cm dish.
7. In 5–7 days, when the culture is confluent, the flask can be subcultured 1:3 using 1 mg/mL collagenase IV or TrypLE Select. A 1X PBS wash is recommended to remove any proteins that may inhibit trypsin activity. The incubation time necessary will be determined by the dissociation agent used.
8. Inactivate the enzyme with an equal volume of complete growth medium and pipet the cells to break up the colonies. The use of a trypsin-like enzyme will allow dissociation of colonies to single cell.
9. Spin the tube at 270 x g for 5 minutes in order to pellet the cells.
10. Resuspend the cells and add to flasks for further expansion.

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)



**UK & Rest of the World**  
184 Park Drive, Milton Park  
Abingdon OX14 4SE, UK  
T: +44 (0)1235 828 200  
F: +44 (0) 1235 820 482



**North America**  
1035 Cambridge Street,  
Cambridge, MA 02141  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218



**Germany**  
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880



**Switzerland**  
Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
T: +41(0) 91 604 55 22  
F: +41(0) 91 605 17 85