


amsbio

CELLBANKER® CRYOPRESERVATION



Survival rate of
all cell types*

CELLBANKER® is a series of easy-to-use cell freezing media offering superior protection against cell stress during freeze/thaw cycles, allowing successful cryopreservation of all mammalian cell types regardless of their sensitivity.

Available in several formulations to suit your research needs, the CELLBANKER® cryopreservation range delivers high rates of cell viability, and reduces post-thaw necrosis and apoptosis, improving cell function and recovery.

CELLBANKER® SERIES

Discover the advantages of CELLBANKER®:

- Consistently high cell viability post thaw (>90%)
- No programmed freezer or liquid nitrogen required
- Long shelf life
- Tested on many cell types
- Freeze directly in -80°C freezer
- Allows viable long term cell storage > 10 years at -80°C or -196°C
- Ready-to-use formulas avoid mixing of DMSO and filtration

CELLBANKER® 1



CELLBANKER® 1 is a ready-to-use serum containing cell cryopreservation medium designed for a broad spectrum of mammalian cells. This unique formulation allows for stable cryopreservation and high viability after freeze-thaw cycles, even for sensitive cells.


CELLBANKER® 2 is a serum-free cell cryopreservation medium optimised for serum-free cultured cells and peptide/protein expressing cells. CELLBANKER® 2 offers a low cost-animal free cryopreservation media for all applications where risks of contamination must be avoided.

CELLBANKER® 2



Typical cell recovery rates for CELLBANKER® 1 and 2

Cell Type	Preservation Period (year)	Viability of cells (%) at -80°C
Fibroblasts	5	90
Caco-2	5	90
K562	5	90
Jurkat	10	90



Viability Comparative Trial (-80°C, 1 Year Storage)

Comparative Analysis of CELLBANKER® 1 vs. conventional cell cryopreservation method.

The unique formulation allows for stable cryopreservation and high viability after freeze-thaw procedure, even for sensitive cell lines.

STEM-CELLBANKER®
GMP grade

Available DMSO Free




STEM-CELLBANKER® is a chemically defined freezing media optimized for stem cells and iPS cells storage, as well as fragile primary cells. Furthermore, recent data supports its ability to cryopreserve organoids and tissues to allow the recovery of viable cells.

STEM-CELLBANKER® GMP grade is manufactured in compliance with JPN, EU, US, and PIC/S GMP guidelines. Free from animal derived components, and contains only chemically defined USP, EP and JP grade ingredients.

A Drug Master File for STEM-CELLBANKER® has been registered with the Food and Drug Administration (FDA) (MF# 15785).

Available in both DMSO containing and DMSO-Free formulations, STEM-CELLBANKER® is an optimal freezing solution for basic research and clinical application of cell therapy products.



STEM-CELLBANKER® vs Conventional freezing

Cell viability assessment: trypan Blue & calcein-esterase based live-dead assay.

Result: Significantly higher cell viability was observed while cell pluripotency, normal karyotype and proliferation ability were maintained in HS293, HS306 and ChiPSA cell lines.

Study by Department of Clinical Science, Intervention and Technology, Karolinska Institute.

STEM-CELLBANKER®
EX GMP grade


New product



STEM-CELLBANKER® EX is a serum free and xeno free, sterile DMSO-containing cryopreservation solution which consists of ingredients already approved in the manufacture of cell therapy excipients.

STEM-CELLBANKER® EX is optimized for cryopreserving human adipocyte-derived mesenchymal / stem cells and human hematopoietic stem cells.

Cryopreservation performance test using Adipose-Derived Stromal/Stem Cells (ADSC) in STEM-CELLBANKER® EX



(A) Cell Viability immediately after thawing cryopreserved cells.

(B) Cell numbers were analyzed after thawing and overnight incubation.

(C) Cell proliferation was evaluated after thawing and culturing for 72 hours.

In summary: High cell viability and maintenance of cell function of STEM-CELLBANKER® EX cryopreserved cells after thawing

HSC-BANKER®



HSC-BANKER® is an optimized GMP grade cryopreservation medium for hematopoietic stem cells.

With a simple formulation HSC-BANKER® also offers the potential to be used as excipient in cell therapy delivery.

Studies show that the HSC-BANKER® is at least equivalent to conventional protocols using DMSO and DEXTRAN.

Featured References

- Nature Methods 2017 May 14:743-751
Fused cerebral organoids model interactions between brain regions.
- Stem Cell Research & Therapy 2016 Jan 12;7:8
Biological impact of xeno-free chemically defined cryopreservation medium on amniotic epithelial cells.
- Stem Cell Research & Therapy 2015 Sep 6:162
Transplantation of cultured dental pulp stem cells into the skeletal muscles ameliorated diabetic polyneuropathy: therapeutic plausibility of freshly isolated and cryopreserved dental pulp stem cells
- Sci Rep 2014 Jan 8;4:3594
A novel efficient feeder-free culture system for the derivation of human induced pluripotent stem cells.
- Nat Protoc. 2014 Oct;9(10):2354-68
Monolayer culturing and cloning of human pluripotent stem cells on laminin-521-based matrices under xeno-free and chemically defined conditions.


Description	Pack Size	Catalogue No.
CELLBANKER® 1 - Serum containing	20 ml	11889
CELLBANKER® 1 - Serum containing	100 ml	11888
CELLBANKER® 2 - Serum free	20 ml	11892
CELLBANKER® 2 - Serum free	100 ml	11891
STEM-CELLBANKER® - GMP	20 ml	11897
STEM-CELLBANKER® - GMP	100 ml	11890
STEM-CELLBANKER® - GMP - DMSO - Free	20 ml	11897F
STEM-CELLBANKER® - GMP - DMSO - Free	100 ml	11890F
STEM-CELLBANKER® EX - GMP	100ml	11890EX
HSC-BANKER® - GMP	15ml	11900
CELLOTION® Cell Wash Solution	100ml	11898

*For full list of cell types tested please contact AMSBIO.


AMSBIO is the global source for CELLBANKER® reagents. CELLBANKER® is a registered trade mark of and manufactured by NIPPON ZENYAKU KOGYO.

AMSBIO | www.amsbio.com | info@amsbio.com

V1

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SL, U.K.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945 5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Via Lisano 3
(CP.683)
6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85